**NEGATOSCOPIO**

**NEGATOSCOPE**

**NEGATOSKOPE**

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**NEGATOSCOPIO**

# **COMPACT 15/30 - 230V**

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*Dear Doctor,*

*Thank you for the faith that you have shown us by acquiring our negatoscope **COMPACT**.*

*We want you to know that our firm is at your complete disposition to furnish information and to answer any questions concerning the operation and use of this equipment.*

*We remind you that for the correct usage of this machine it is necessary to first read this manual. The **COMPACT** corresponds to all of the relative safety regulations and presents no danger to the operator when used according to the instructions.*

*While wishing you success in your work, we remind you that the reproduction of this manual is forbidden and that the technical components could be changed without warning due to our continued technical research.*

**EURONDA S.p.A.**

## • GENERAL INDEX •

• <b>INTRODUCTION</b>	
GENERAL NOTES ON DELIVERY	P. 3
GENERAL INSTRUCTIONS	P. 3
MEANING OF SYMBOLS	P. 4
WARRANTY INFORMATION	P. 4
• <b>PART 1</b>	
MACHINE AND PACKING OVERALL DIMENSIONS	P. 5
• <b>PART 2</b>	
INSTALLATION	P. 6
INSTRUCTIONS FOR WALL-MOUNTING	P. 6
SAFETY	P. 7
ELECTRICAL CONNECTION	P. 8
• <b>PART 3</b>	
DESCRIPTION	P. 9
TECHNICAL AND OPERATING CHARACTERISTICS	P. 9
• <b>PART 4</b>	
OPERATING INSTRUCTIONS	P.10
• <b>PART 5</b>	
ORDINARY MAINTENANCE	P. 11
EXTRAORDINARY MAINTENANCE AND DIRECTIVES	P. 11
REPLACEMENT OF COMPONENTS	P. 11
WARNINGS FOR THE USER	P. 12
PROBLEM SOLVING	P. 13
• <b>PART 6</b>	
WIRING DIAGRAMS	P. 14
- Wiring Diagram COMPACT 15	P. 14
- Wiring Diagram COMPACT 30	P. 14
NEGATOSCOPE COMPONENTS	
- <b>COMPACT 15</b>	
- Top view/Bottom view	P. 15
- Inside view	P. 16
- <b>COMPACT 30</b>	
- Top view/Bottom view	P. 17
- Inside view	P. 18
SCRAPPING INSTRUCTIONS	P. 19

## GENERAL NOTES ON DELIVERY

Upon receipt of the equipment, please control the condition of the package and keep it for possible shipments.

Open the package and check that:

- the equipment corresponds to the specific requirements
- there is no evident damage.

In case of damage or missing parts, inform immediately either the shipping agent, the dental depot or Euronda S.p.A.

All of the dimensions listed in this manual are not binding.

The sketches and any other consigned documents regarding this machine are property of Euronda S.p.A., which reserves all rights. They can not be made available to any third parties.

It is also forbidden to reproduce in whole or in part either the manual or the illustrations.

## GENERAL INSTRUCTIONS

Included in this manual are instructions for:

- correct installation
- the safe and efficient use of the apparatus
- routine maintenance.

The COMPACT must be used in compliance with the instructions given in this manual and never for other purposes.

The user is legally responsible for the relative performance of the installation, operation and use of this apparatus.

If the apparatus is not operated correctly or maintained adequately, the manufacturer is not responsible for the eventual breakage, damage, or faulty operation of the equipment.

**Use of the machine:** this manual is an integral part of the ultrasound apparatus and it is necessary to keep it near to the machine in order to ensure an easy reference.

The machine is for professional use only.

Use the machine only for the ultrasound functions outlined in this manual.

**COMPACT - 230V NEGATOSCOPE: instrument designed and developed for the examination of endoral and/or panoramic radiographies.**

## MEANING OF SYMBOLS

The instrument label shows a number of symbols, as prescribed by the UNI EN 980 standard. Their respective meaning is specified below.

SYMBOL	DESCRIPTION
<b>SN</b>	<b>"SERIAL NUMBER"</b> This symbol must be accompanied by the manufacturer's serial number. The serial number must be adjacent to the symbol.
	<b>"DATE OF MANUFACTURE"</b> This symbol must be accompanied by the year. The year must be written as a four-figure number.
	<b>"ATTENTION, SEE OPERATING INSTRUCTIONS"</b>

Moreover, in the present manual, various symbols can be seen at the side of descriptions, notes, etc. The purpose of these illustrations is to draw the reader's attention to a particular note or explanation. The meaning is specified below.

SYMBOL	DESCRIPTION
	<b>PARTICULARLY IMPORTANT INFORMATION FOR REASONS OF SAFETY</b> This symbol is intended to draw the reader's attention to specific fundamental principles for the operator's safety.

## GUARANTEE INFORMATION

Euronda guarantees the quality of its equipment, if used in accordance with the instructions supplied in this manual, and according to the conditions printed on the GUARANTEE CERTIFICATE.

**WARNING:** the CUSTOMER USER must fill in all parts of the DETACHABLE SLIP of the guarantee certificate and send it to EURONDA S.p.A.


The guarantee is effective from the delivery date of the unit to the customer, which is confirmed by the return of the guarantee slip correctly filled in and signed. In the case of dispute, the date indicated on the purchase invoice, which shows the serial number of the unit, will be considered valid.

**NOTE:** Please retain the original packaging and use it each time the unit needs to be transported. The use of different packaging could damage the product during transportation.

**The unit will not be replaced nor will the warranty be extended in the event of misuse. Repairs under warranty will be effected according to the unquestionable decision of Euronda S.p.A. and will not include packaging, transport or labour costs.**

**PART 1**  
**MACHINE AND PACKING OVERALL DIMENSIONS**

English

**COMPACT 15**

**COMPACT 30**

<b>Dimensions (mm):</b> <b>W= width; H= height; D= depth</b>	<b>COMPACT 15</b>	<b>COMPACT 30</b>
MACHINE (W x H x D)	420 x 230 x 71	495 x 380 x 71
PACKING (W x H x D)	540 x 300 x 150	615 x 450 x 150

## PART 2


### INSTALLATION

**The COMPACT negatoscope has been tested and inspected in the factory. This instrument does not create electromagnetic interference with other devices.**

Remove the packing and install the instrument observing the following instructions:

1. Position the negatoscope on a flat surface or on the wall.
2. Do not install the machine close to a sink or similar to avoid contact with water and possible electrical short circuits.
3. Keep the machine far from heat sources
4. Ensure that the electrical cable is not folded up and freely reaches the electrical socket.
5. Connect the supply cable socket to the negatoscope recessed plug, and then connect the plug to the mains: the instrument is now ready for use.

### INSTRUCTION FOR WALL MOUNTING


1. Drill two holes in the wall, keeping "A" distance between centres.

**COMPACT 15: A = 291 mm**

**COMPACT 30: A = 200 mm**

Use a suitable dowel, depending on the type of wall (cement, brick, etc.).

2. Rest the negatoscope against the wall in order to insert the two screw heads into the appropriate ∅ 7 holes in the rear hinges.
3. Gently slide the negatoscope downwards, until it locks into position.

## SAFETY

**Before using your COMPACT negatoscope, carefully read these safety warnings. Failure to observe them can cause accidents or damage the apparatus.**

**KEEP THE AREA AROUND THE NEGATOSCOPE CLEAN AND DRY.**

**NEVER DISASSEMBLE THE COMPACT® NEGATOSCOPE.  
THE HIGH INTERNAL VOLTAGE IS DANGEROUS**

If it is not possible to disconnect the electricity or if it is too far or not visible to the person who is making the repairs, attach the sign "Work in Process" on the switch after having turned it "OFF".

**CLEAN THE NEGATOSCOPE WITH A DAMP CLOTH, AFTER ENSURING THAT THE CONNECTION CABLE HAS BEEN UNPLUGGED (before reusing the machine make sure that the machine is completely dry).  
DO NOT USE SOLVENTS IN THE LABEL.**

**MAKE SURE THAT THE WIRING HAS A GROUND CABLE WHICH CONFORMS TO ALL OF THE SAFETY NORMS OF THE COUNTRY OF INSTALLATION.**

**NEVER REMOVE THE LABEL ON THE MACHINE**

**ENSURE THAT THE MACHINE IS CONNECTED TO THE CORRECT ELECTRICAL VOLTAGE**

## USE EXCLUSIVELY THE ORIGINAL REFILLS

If the negatoscope doesn't WORK properly refer to the selection in the manual "Solution to Problems" for the possible cause. If you need subsequent information or the machine have to be repaired contact the specialised personnel of the supplier or the technical office of Euronda S.p.A.

## WARNING

- The negatoscope conforms to the electrical safety requirements set out by the Standards Institute, and is provided with a bipolar plug which ensures complete earthing.
- As a fundamental safety requirement, make sure that the electric installation is efficiently earthed and the capacity of the installation and electric current sockets is adequate for the apparatus power shown on the rating plate. Have qualified personnel carry out the above checks.

**BY NOT OBSERVING THE REQUIREMENTS LISTED EURONDA S.p.A. FORFEITS ALL RESPONSABILITY FOR ITS MISUSE.**

## ELECTRICAL CONNECTION

Check that the electrical current indicated on the tag on the rear of the machine, corresponds to the voltage of the place of installation.

The machine must be connected to an electrical system with adequate voltage which corresponds to all standards laws of the country of installation.

The **earth** terminal **protecting** the instrument can be connected to the system earth wire.


In fact, there is a fixed and permanently fitted M4 through screw, which is internally connected to the negatoscope protection earthing, and can be externally connected to the system earth wire.

To connect the system earth wire:

- 1) unscrew the knurled nut positioned on the bottom part of the COMPACT negatoscope, and take out the toothed washer;
- 2) insert the system earth wire terminal;
- 3) insert the toothed washer and retighten the knurled nut.

Knurled nut

**ATTENTION:** Don't fold or twist the electrical cable and don't lay any objects on it.

## PART 3

### DESCRIPTION

**COMPACT:** universal table or wall-mounted negatoscope in ABS, suitable for all the requirements of a dental consulting room. It is provided with two or four white-light fluorescent tubes and a radiography holder for the examination of standard radiographies. Its compact design and solid structure ensure excellent fluorescent luminosity with a continuous luminous flux.


### TECHNICAL AND OPERATING CHARACTERISTICS

CARATTERISTICHE	COMPACT 15	COMPACT 30
Power supply voltage	230 Volt	230 Volt
Power supply frequency	50 Hz	50 Hz
Number of phases	1	1
Absorbed power	44 W	95 W
Insulation class	I	I
No. of ballasts	1	2
Type	30 / 32 W - 230 V	30 / 32 W - 230 V
No. of neon lights	2	4
Type of neon	F14W/154-T8	F15W/154-T8
Type of light	6500 K	6500 K
Socket	G13	G13
Dimensions (mm)	Ø26 L = 360	Ø26 L = 438
Switch	Bipolar	Bipolar
Fuses	N° 2 - 5 x 20 - T 0,8 A	N° 2 - 5 x 20 - T 2 A
Weight	2,6 kg	4,6 kg
Useful surface of the viewer (mm)	340 x 145	416 x 300
Colour	RAL 9002	RAL 9002


## PART 4

### OPERATING INSTRUCTIONS

1. Make sure that the negatoscope has been electrically connected.  
**Press the lateral POWER switch placed on the side of the instrument.**  
The negatoscope screen will reflect the white light required.

2. Insert the radiography between the radiography holder and the Plexiglas viewer, as shown in the illustration below.

3. The radiography can now be examined properly.

## PART 5

### ORDINARY MAINTENANCE

- The negatoscope does not require any ordinary maintenance.
  - Disconnect the main supply cable from the plug (mains voltage).
  - Use a slightly damp cloth to wash the outside of the instrument
- Attention on the opal viewer use only water or neutral detergents; do not use alcohol, disinfectants, etc...

**Ensure that the negatoscope is dry before the switch on.**

### EXTRAORDINARY MAINTENANCE AND DIRECTIVES

- **ALL THE MAINTENANCE OPERATIONS, INCLUDING SUPPLY CABLE REPLACEMENT, MUST BE CARRIED OUT BY COMPETENT PERSONNEL.**
- Each extraordinary maintenance operation must be carried out as specified in PART 2.

### REPLACEMENT OF COMPONENTS (fluorescent lamp – starting switch - Fuses)

#### **ATTENTION:**


**BEFORE EVERY MAINTENANCE INTERVENTION TURN OFF THE VOLTAGE AND ALLOW A SPECIALISED TECHNICIAN TO PERFORM THE CHECK**

#### **DANGER:**

**HIGH INTERNAL VOLTAGE**

To replace faulty components, open the instrument as follows: remove the fixing screws and lift up the top guard.

The various components will now be accessible.

#### **To replace the fluorescent lamp tubes:**

1. rotate the faulty tube on its axis by approximately a quarter turn (a click should be heard);
2. take the tube out of the lamp socket;
3. insert a new tube in the lamp socket;
4. rotate the tube as specified in point 1.

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#### **To replace the fuses:**

1. take out the fuse carrier;
2. take out the faulty fuse from its seat;
3. insert a new fuse;
4. insert the fuse carrier in its appropriate housing.

- **You are always recommended to use spare parts having the same technical characteristics. When requesting spare parts, always give the code specified in the list of components.**

### WARNINGS FOR THE USER

- The installation and replacement of lamp sockets, starting-switch holders, starting switches, tubes and fuses must be carried out by qualified personnel.
- Never leave a lamp socket without its lamp, as some of the parts may be live; **therefore disconnect the electricity before replacing the lamp.**
- The fluorescent lamps which are exhausted or about to be exhausted must be promptly replaced, since their malfunction can cause overheating in the lamp socket or starting-switch holder.
- In the case where the lamp sockets show obvious deterioration, have them replaced and checked to find out the cause of the faults which may have produce the damage.

## PROBLEM SOLVING

PROBLEMS	CAUSES	SOLUTIONS
The instrument does not light up	No voltage	Wait and/or restore voltage
	Protection fuses	Check any short circuits or faults, then replace with fuses having the same value
	Faulty ballast	Replace the ballast
	Faulty fluorescent lamp	Replace the lamp
Blinking light	The gas flux in the tube is not uniform during the first few hours of operation	Leave the instrument on for 7-8 hours
	Line voltage too low	Check the voltage in the Surgery
	Defective fluorescent lamp	Replace it with one having the same characteristics
The fluorescent lamp keeps switching on and off	Faulty starting switch	Replace the starting switch
Areas of shadow are left on the luminous surface	Exhausted fluorescent lamp	Replace the fluorescent lamp with another having the same characteristics


**ATTENTION: BEFORE EVERY MAINTENANCE INTERVENTION  
TURN OFF THE VOLTAGE ALLOW A SPECIALISED  
TECHNICIAN TO PERFORM THE CHECK**


**DANGER: HIGH INTERNAL VOLTAGE**

# **PART 6** **ELECTRICAL DIAGRAMS**

## **Negatoscope model COMPACT 15 - 230 V**

## **Negatoscope model COMPACT 30 - 230 V**


## NEGATOSCOPE COMPONENTS COMPACT 15

Top view

Bottom view

# Inside view

English

**COMPACT 30****Top view****Bottom view**

# Inside view

## SCRAPPING INSTRUCTIONS

The COMPACT® negatoscope has been manufactured using iron materials, electronic components and plastic materials. No particular precautions are required to throw the unit away. Do not dump the unit in an unguarded place; consign it to a waste disposal agency for scrapping. Always refer to the specific laws in the country of use when scrapping equipment.

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Copyright © 1999 Euronda S.p.A.